

In pharmaceutical manufacturing, precision and hygiene are paramount. Valves help achieve these by regulating the flow of ingredients, maintaining sterile conditions, and preventing contamination. They are used in processes such as mixing, batching, transferring, and filling, where precise control over the flow and pressure of fluids is essential. The reliability and performance of valves directly impact the quality of the final product and the efficiency of the manufacturing process.

## **Company Background**

Since 1969, AVCON has been producing Solenoid Valves, Angle Seated Valves, Control Valves, Electric and Pneumatic Actuators, Ball and Butterfly Valves. However, we started manufacturing products for the pharmaceutical sector in 2001.

When the pharmaceutical business offered creative valve makers the chance to domestically design and develop diaphragm valves to regulate SIP, CIP, and purified water system, we began producing hygienic diaphragm valves in 2001. With highly dependable domestically produced diaphragm valves that satisfied industry criteria for material, drainability, particle minimization, surface quality, and the accuracy requirements of low offset, AVCON was the first Indian company to enter the domestic market. Since then, AVCON Valves have been the most often specified for mixing, batching, transferring, filling, sterile filtration, purified water, cosmetics, and process water.

## **Our Synopsis**

AVCON was established to support ongoing innovation and expand the range of products we provide. Since 1969, we have successfully installed over 2,50000 industrial process valves and automation for the steel, cement, textile, chemical, and food and over 1,00000 valve installations in pharmaceuticals industry segments.



## **Quality Objectives**

Our ultimate objective at AVCON is to supply valves with lower operating and maintenance costs while increasing our customers' asset utilization and product throughput.

Invest in our Customers, goods, and services while implementing ideas and creating solutions for our clients at a purposeful pace.

## Industries that we serve

Avcon Controls Pvt. Ltd. has walked through several Years, gaining valuable experience, offering Avcon equipments in the Products Manufactured by the OEMs, in various industries such as Food, Beverages, Pharmaceutical, Biotech, Cosmetic, Personal Care, Dairy, Mechanical, Chemical, Food Processing, Water Treatment & Purification, Wine Production & for all types of clean or sterile processing applications.









# Manual Diaphragm Valve 6227M / 6227F Series

#### **Design & Manufacturing Standards**

Design Standard ASME BPE
Testing Standard EN12266-1

Certification Standard EN-10204:2.1,2.2,3.1
Surface Finish ASME BPE (SF4)
Diaphragm Compliance USP CLASS VI Biological
Reactivity Test In Vitro and In

Vivo US FDA CFR177.1550 (PTFE) US FDA CFR177.2600

(EPDM)

#### **Benefits**

Precise Flow Control Bubble-Tight Sealing Resistance to Corrosion

**Excellent Containment Properties** 

Biocompatibility

Smooth Flow Characteristics

Low Pressure Loss

Elastomeric Diaphragms for Enhanced Performance

Easy Automation

Reduced Water Hammer

Long Service Life

SIP and CIP compatibility

#### **Applications**

Mixing

Batching

Transferring

Filling

Sterile filtration

CIP & SIP

Cosmetics

Purified Water & Process Water

#### **Operating Conditions**

Service Mixing, Batching, Transferring, and

Filling, Purified Water & Process Water

Operating Temp Ambient/85°C

Design Temp 150 ° C
Operating Pr. 10 Kg/cm²(g)
Design Pressure 12 Kg/cm² (g)/FV

Vapour Yes

Viscosity Cp As of service

#### **Material Of Construction**

Body CF3M (Investment casted and

machined) SS-316L (Forged &

Machined)

Stem Assembly AISI 304

Diaphragm Material EPDM Or PTFE with EPDM backup

Bonnet AISI 304

Hand wheel 30% graphited Nylon (PA-6)

Packing/Sealing PTFE Studs/Nuts SS 304





#### **End Connections**

Inlet Butt weld/TC end/Flanged (Optional)
Outlet Butt weld/TC end/Flanged (Optional)

#### **Finish**

 $\begin{array}{ll} \text{Internal} & \text{Electro polished to } < 0.38 \mu \text{m Ra} \\ \text{External} & \text{Mirror Polished to } 180 \, \text{Grit} \end{array}$ 

#### **Performance Characteristics**

 $\begin{array}{ll} \mbox{Hydro Test Pressure} & \mbox{Body 15.0 Kg/cm}^2(g) \\ \mbox{Hydro Test Pressure} & \mbox{Seat 10Kg/cm}^2(g) \end{array}$ 




Stroke Limiter SS 304 with Plastic cap

## **Manual Diaphragm Valve**

## 6227M Series

## 6227M - Casted SPG






CAT NO.	SIZE	ØA	ØB	С	F	T
BS6227M086	1/4"	Ø13.5	Ø40	68.5	70	1.65
BS6227M166	1/2"	Ø12.7	Ø58	76	112	1.65
BS6227M216	3/4"	Ø19.05	Ø80	104	117	1.65
BS6227M266	1"	Ø25.4	Ø80	104	127	1.65
BS6227M416	1-1/2"	Ø38.1	Ø118	142	146	1.65
BS6227M516	2"	Ø50.8	Ø118	144	159	1.65
BS6227M666	2-1/2"	Ø63.5	Ø118	160	190	1.65
BS6227M816	3"	Ø76.2	Ø200	250	270	1.65
BS6227M1016	4"	Ø101.6	Ø200	291.7	292.1	2.11

## 6227M - Casted TC






CAT NO.	SIZE	ØA	ØВ	С	F	Т
BS6227M167	1/2"	Ø25	Ø58	76	112	2.85
BS6227M217	3/4"	Ø50.5	Ø80	104	117	2.85
BS6227M267	1"	Ø50.5	Ø80	104	127	2.85
BS6227M417	1-1/2"	Ø50.5	Ø118	142	146	2.85
BS6227M517	2"	Ø64	Ø118	144	159	2.85
BS6227M667	2-1/2"	Ø77.4	Ø118	160	190	2.85
BS6227M817	3"	Ø90.9	Ø200	250	270	2.85
BS6227M1017	4"	Ø118.9	Ø200	291.7	292.1	2.85

## 6227F - Forged SPG






CAT NO.	SIZE	ØA	ØB	С	F	Т
BS6227F166	1/2"	Ø12.7	Ø58	75	112	1.65
BS6227F216	3/4"	Ø19.05	Ø80	103	117	1.65
BS6227F266	1"	Ø25.4	Ø80	103	127	1.65
BS6227F266	1"	Ø25.4	Ø80	116	127	1.65
BS6227F416	1-1/2"	Ø38.1	Ø118	136	146	1.65
BS6227F516	2"	Ø50.8	Ø118	144	159	1.65
BS6227F666	2-1/2"	Ø63.5	Ø118	160	190	1.65

## 6227F - Forged TC





CAT NO.	SIZE	ØA	ØE	ØВ	С	F	M
BS6227F167	1/2"	Ø25	Ø9.4	Ø40	68.5	90	2.85
BS6227F167	1/2"	Ø25	Ø9.4	Ø58	75	112	2.85
BS6227F217	3/4"	Ø50.5	Ø15.75	Ø80	103	117	2.85
BS6227F267	1"	Ø50.5	Ø22.1	Ø80	103	127	2.85
BS6227F467	1-1/2"	Ø50.5	Ø34.8	Ø118	136	146	2.85
BS6227F517	2"	Ø63.9	Ø47.5	Ø118	144	159	2.85
BS6227F667	2-1/2"	Ø77.4	Ø60.2	Ø118	160	190	2.85

ā avcon<sup>®</sup>

## **Pneumatic Diaphragm Valve**

### **6827M/6827F Series**

#### **Design & Manufacturing Standards**

Design Standard ASME BPE
Testing Standard EN12266-1

Certification Standard EN-10204: 2.1,2.2,3.1 Surface Finish ASME BPE (SF4)

Diaphragm Compliance USP CLASS VI Biological Reactivity

Test In Vitro and In Vivo US FDA CFR177.1550 (PTFE) US FDA CFR177.2600 (EPDM)

#### **Benefits**

Precise Flow Control
Bubble-Tight Sealing
Resistance to Corrosion

**Excellent Containment Properties** 

Biocompatibility

**Smooth Flow Characteristics** 

Low Pressure Loss

Elastomeric Diaphragms for Enhanced Performance

Easy Automation Reduced Water Hammer Long Service Life

SIP and CIP compatibility

#### **Applications**

Mixing Batching Transferring

Filling

Sterile filtration CIP & SIP

Cosmetics
Purified Water & Process Water

#### **Material of Construction**

Body CF3M (Investment casted and machined)

 $SS-316L \, (Forged \, \& \, Machined)$ 

Stem Assembly AISI 304

Diaphragm Material EPDM Or PTFE with EPDM backup

Bonnet AISI 304

Hand wheel 30% graphited Nylon (PA-6)

Packing/Sealing PTFE Studs/Nuts SS 304

#### **Pneumatic Actuator**

Type Linear

Construction Piston actuated

Material of construction Plastic (Gr Pa6) & SS 304

Input

Air Signal Supply Bar (g) 4 to 6 BAR Pneumatic Connections 1/4" BSP (F)

Fail Safe Position Normally Closed / Open

Local Position Indicator Impact Resistant Thermoplastic Dome

#### **End Connections**

 $\label{thm:continuous} In let \quad Spigot/Triclover\,end\,as\,per\,BS\,4825/ASME\,BPE/Flanged\,(Optional)\\ Outlet \quad Spigot/Triclover\,end\,as\,per\,BS\,4825/ASME\,BPE/Flanged\,(Optional)\\$ 

#### **Finish**

 $\begin{array}{ll} \text{Internal} & \text{Electro polished to } < 0.38 \mu\text{m Ra} \\ \text{External} & \text{Mirror Polished to } 180 \, \text{Grit} \end{array}$ 

#### **Performance Characteristics**

Hydro Test Pressure Body 15.0 Kg/ cm²(g) Hydro Test Pressure Seat 10Kg/cm²(g)



### **Accessories**

#### Solenoid Valves (optional)

Type 3/2 way, Normally close energies to open Coil voltage 230AC/24 DC/suitable for safe area

Pressure Bar (g) 10 BAR Max Material PA6/Brass/SS-304 Air Connection 1/4" BSPT (F)

Cable Entry Plug in connector PG-9

Area Classification Safe area, Weatherproof suitable for min. IP 65

#### **Proximity Sensor** (optional)

Type of Switch Inductive Proximity or Limit switch

Area classification IP 65 No. of Switches 2 Nos.

Cable Entry Supply Voltage 24VDC

#### Stroke Limiter (optional)

Stroke Limiter SS 304 with Plastic cap


#### **Electronic Valve Positioner** (Optional)

MOC Plastic (Gr PA6-GF30) (Standard) SS-304 (Optional)

Signal Input 4-20mA
Power Supply 24VDC
Protection Class IP-66
Ambient Temperature 0~70°C
Explosion-Proof Class Ex-nA IICT4
Signal air Pressure 4-6 Bar


# Pneumatic Diaphragm Valve 6827M (Investment Casted)

## 6827M - Casted SPG & TC Pneumatic Diaphragm Seated Valves (Pa6).




	Valve	Model No	Actuator (normally	Si	ze	Flow Kv	Pressure	Drain Angle	Weig ht		[	Dimen	sions		
	TC	SPG	close)	inch	mm	МЗНг	Kg/cm 2	u	Kgs	ØΑ	F	Р	Q	Т	ØВ
	6827M167/TE/TC	6827M166/TE/SPG	K-52	1/2"	15	3.5	10.0	20°	0.7	12.7	112	124	63.5	1.65	25
	6827M217/TE/TC	6827M216/TE/SPG	K-62	3/4"	20	5.0	10.0	21.5°	1.2	19.05	117	150	80	1.65	50.5
	6827M267/TE/TC	6827M266/TE/SPG	K-62	1"	25	12.5	10.0	21.5°	1.3	25.4	127	150	80	1.65	50.5
	6827M417/TE/TC	6827M416/TE/SPG	K-102	1-1/2"	40	36.0	10.0	25°	5	38.1	146	240	131	1.65	50.5
	6827M517/TE/TC	6827M516/TE/SPG	K-122	2"	50	42.0	10.0	30°	7.3	50.4	159	282	161	1.65	63.9
1	6827M667/TE/TC	6827M666/TE/SPG	K-123	2-1/2"	65	53.0	6.0	31.5°	8.1	63.5	190	300	161	1.65	77.4


## 6827M - Casted SPG & TC Pneumatic Diaphragm Seated Valves (SS).



Valve	Model No	Actuator (normally	Si	ze	Flow Kv	Pressure	Drain Angle	Weig ht			Dimer	sions	6	
TC	SPG	close)	inch	mm	МЗНr	Kg/cm 2	u	Kgs	ØA	F	Р	Q	Т	ØВ
N.A.	6827M086/TE/SPG	M-42	1/4"	8	3.5	10	35°	0.7	12.7	90	100	47	1.65	N.A.
6827M167/TE/TC	6827M166/TE/SPG	M-52	1/2"	15	3.5	10	50°	1.3	12.7	112	124	59	1.65	25
6827M217/TE/TC	6827M216/TE/SPG	M-62	3/4"	20	5	10	45°	2.5	19.05	117	171.5	74	1.65	50.5
6827M267/TE/TC	6827M266/TE/SPG	M-62	1"	25	12.5	10	40°	2.6	25.4	127	171.5	74	1.65	50.5
6827M417/TE/TC	6827M416/TE/SPG	M-102	1-1/2"	40	36	10	35°	10.2	38.1	146	273	120	1.65	50.5
6827M517/TE/TC	6827M516/TE/SPG	M-122	2"	50	42	10	35°	14.0	50.4	159	300	144	1.65	63.9
6827M667/TE/TC	6827M666/TE/SPG	M-123	2-1/2"	65	53	06	30°	15.0	63.5	190	312	144	1.65	77.4
6827M817/TE/TC	6827M816/TE/SPG	CY-150	3"	80	62	06	31°	17.9	90.9	270	375	180	1.65	90.9
6827M1007/TE/TC	6827M1006/TE/SPG	CY-150	4"	100	72	06	31°	19.9	101.6	292	425	190	2.11	118.9


### **Pneu. Valve with valve Positioner**






Valve Model No	Actuator (Normally	Siz	е	Dimensions				
PTFE/EPDM	Close)	inch	mm	В	С	F		
6827M166/167	K-52	1/2"	15	63.5	199	112		
6827M216/217	K-62	3/4"	20	80	225	117		
6827M266/267	K-62	1"	25	80	225	127		
6827M416/417	K-102	1 1/2""	40	131	315	146		
6827M516/517	K-122	2"	50	161	357	159		
6827M666/667	K-123	2 1/2"	65	161	375	190		
6827M816/817	CY-150	3"	80	180	450	270		
6827M1016/1017	CY-150	4"	100	213	517	292		

## Pneu. Valve with Feedback Unit






Valve Model No	Actuator (Normally	Size	9	Dimensions				
PTFE/EPDM	Close)	inch	mm	В	С	F		
6827M166/167	K-52	1/2"	15	63.5	199	112		
6827M216/217	K-62	3/4"	20	80	225	117		
6827M266/267	K-62	1"	25	80	225	127		
6827M416/417	K-102	1 1/2""	40	131	315	146		
6827M516/517	K-122	2"	50	161	357	159		
6827M666/667	K-123	2 1/2"	65	161	375	190		
6827M816/817	CY-150	3"	80	180	450	270		
6827M1016/1017	CY-150	4"	100	213	517	292		

'avcon<sup>®</sup>


# Pneumatic Diaphragm Valve 6827F (Forged)

## 6827F - Forged SPG & TC Pneumatic Diaphragm Seated Valves (Pa6).




	Valve	Model No	Actuator (normally	Si	ze	Flow Kv	Pressure	Drain Angle	Weig ht			Dimen	sions		
	TC	SPG	close)	inch	mm	МЗНr	Kg/cm 2	u	Kgs	ØΑ	F	Р	Q	Т	ØВ
	6827F167/TE/TC	6827F166/TE/SPG	K-52	1/2"	15	3.5	10.0	20°	0.7	12.7	112	124	63.5	1.65	25
	6827F217/TE/TC	6827F216/TE/SPG	K-62	3/4"	20	5.0	10.0	21.5°	1.2	19.05	117	150	80	1.65	50.5
4	6827F267/TE/TC	6827F266/TE/SPG	K-62	1"	25	12.5	10.0	21.5°	1.3	25.4	127	150	80	1.65	50.5
1	6827F417/TE/TC	6827F416/TE/SPG	K-102	1-1/2"	40	36.0	10.0	25°	5	38.1	146	240	131	1.65	50.5
•	6827F517/TE/TC	6827F516/TE/SPG	K-122	2"	50	42.0	10.0	30°	7.3	50.4	159	282	161	1.65	63.9
	6827F667/TE/TC	6827F666/TE/SPG	K-123	2-1/2"	65	53.0	6.0	31.5°	8.1	63.5	190	300	161	1.65	77.4


## 6827F - Forged SPG & TC Pneumatic Diaphragm Seated Valves (SS).



	Valve	Model No	Actuator (normally	Si	ze	Flow Kv	Pressure	Drain Angle	Weig ht			Dimer	sions	3	
	TC	SPG	close)	inch	mm	МЗНг	Kg/cm 2	u	Kgs	ØA	F	Р	Q	Т	ØВ
	6827F087/TE/TC	6827F086/TE/SPG	M-42	1/4"	8	3.5	10	35°	0.7	12.7	90	100	47	1.65	25
	6827F167/TE/TC	6827F166/TE/SPG	M-52	1/2"	15	3.5	10	50°	1.3	12.7	112	124	59	1.65	25
	6827F217/TE/TC	6827F216/TE/SPG	M-62	3/4"	20	5	10	45°	2.5	19.05	117	171.5	74	1.65	50.5
	6827F267/TE/TC	6827F266/TE/SPG	M-62	1"	25	12.5	10	40°	2.6	25.4	127	171.5	74	1.65	50.5
	6827F417/TE/TC	6827F416/TE/SPG	M-102	1-1/2"	40	36	10	35°	10.2	38.1	146	273	120	1.65	50.5
	6827F517/TE/TC	6827F516/TE/SPG	M-122	2"	50	42	10	35°	14.0	50.4	159	300	144	1.65	63.9
•	6827F667/TE/TC	6827F666/TE/SPG	M-123	2-1/2"	65	53	06	30°	15.0	63.5	190	312	144	1.65	77.4
-	6827F817/TE/TC	6827F816/TE/SPG	CY-150	3"	80	62	06	31°	17.9	90.9	270	375	180	1.65	90.9
	6827F1007/TE/TC	6827F1006/TE/SPG	CY-150	4"	100	72	06	31°	19.9	101.6	292	425	190	2.11	118.9


### **Pneu. Valve with valve Positioner**





Valve Model No	Actuator (Normally	Siz	е	Dir	ns	
PTFE/EPDM	Close)	inch	mm	В	С	F
6827F166/167	K-52	1/2"	15	63.5	199	112
6827F216/217	K-62	3/4"	20	80	225	117
6827F266/267	K-62	1"	25	80	225	127
6827F416/417	K-102	1 1/2"	40	131	315	146
6827F516/517	K-122	2"	50	161	357	159
6827F666/667	K-123	2 1/2"	65	161	375	190
6827F816/817	CY-150	3"	80	180	450	270
6827F1016/1017	CY-150	4"	100	213	517	292

## Pneu. Valve with Feedback Unit



Valve Model No	Actuator (Normally	Size	Э	Din	nensio	ns
PTFE/EPDM	Close)	inch	mm	В	С	F
6827F166/167	K-52	1/2"	15	63.5	199	112
6827F216/217	K-62	3/4"	20	80	225	117
6827F266/267	K-62	1"	25	80	225	127
6827F416/417	K-102	1 1/2""	40	131	315	146
6827F516/517	K-122	2"	50	161	357	159
6827F666/667	K-123	2 1/2"	65	161	375	190
6827F816/817	CY-150	3"	80	180	450	270
6827F1016/1017	CY-150	4"	100	213	517	292

#### **Design & Manufacturing Standards**

Design Standard ASME BPE
Testing Standard EN12266-1

Certification Standard EN-10204: 2.1,2.2,3.1 Surface Finish ASME BPE (SF4)

PTFE Compliance US FDA CFR177.1550 (PTFE)

#### **Benefits**

High Flow Rate
Minimal Pressure Loss
Reducing Flow Impedance
High Temperatures
Large Flow Rates
Anti-water Hammer Proof

Used With Contaminated Fluids

High Cycle Life

Optical Position Indication Standard

#### **Applications**

Food and Beverage Pharmaceuticals Chemical Processing Water Treatment HVAC Systems.

Suitable for neutral, Slightly Highly Aggressive and Aggressive Media.

#### **Operating Conditions**

Service Water/Air/Steam/Aggressive

Media

 $\begin{array}{ll} \mbox{Design Temp} & 180\,^{\circ}\mbox{C} \\ \mbox{Operating Pr.} & 10\,\mbox{Kg/cm}^{2}(\mbox{g}) \\ \mbox{Design Pressure} & 12\,\mbox{Kg/cm}^{2}(\mbox{g})/\mbox{FV} \end{array}$ 

Vapour Yes

Viscosity cP As of service

#### **Material of Construction**

Body ASTM 351-CF8/CF8M/CF3M Stem Assembly ASTM 351-CF8/CF8M/CF3M

Seat Material Food Grade PTFE

Studs/Nuts SS 304

#### **Pneumatic Actuator**

Type Linear

Construction Piston actuated (Single

acting or double acting)

Hand Wheel -

 $Material \, of \, construction \qquad Plastic (Gr \, PA \, 6) \, \& \, SS \, 304$ 

Input

Air Signal Supply Bar (g) 4 to 6 BARPneumatic Connections 1/4" BSP (F)

Fail Safe Position Normally Closed / Open

 ${\sf Local \, Position \, Indicator} \qquad {\sf Impact \, Resistant \, Thermoplastic \, Dome}$ 

#### **End Connections**

Inlet Spigot/Triclover (Standard)/
Threaded OR Flanged (Optional)
Outlet Spigot/Triclover (Standard)/

Threaded OR Flanged (Optional)

#### **Finish**

Internal Electro polished to < 0.38 µm

Ra as per ASME BPE (SF4)

External Mirror Polished to 180 Grit

#### **Performance Characteristics**

Hydro Test Pressure Body 15.0 Kg/ cm²(g) Hydro Test Pressure Seat 10Kg/cm²(g)

# Pneumatic Angle Seat Valve 6812M/6212M Series



## **Accessories**

#### Solenoid Valves (optional)

Type 3/2 way, Normally close energise to open Coil voltage 230AC/24 DC/suitable for safe area

Air Supply Bar (g) 4 to 8

Material PA6/Brass/SS-304
Air Connection 1/4" BSPT (F)
Cable Entry Plug in connector PG-9

Area Classification Safe area, Weatherproof suitable for min. IP 65

#### **Proximity Sensor** (optional)

**Proximity Switch** 

Make P&For equivalent

Type Inductive Proximity switches, 2 wire

Area classification Safe area, Weatherproof suitable for min. IP 65 No. of Switches 2 Nos. (1 No. Each for open & close positions)

Cable Entry 1/2" NPT (F) - 2 Nos. (1 No. Plugged)

#### **Electronic Valve Positioner** (Optional)


MOC Plastic (Gr PA6-GF30) (Standard) SS-304 (Optional)

Signal Input 4-20mA
Power Supply 24V DC ± 10%
Ambient Temperature 0-70 ° C
Supply air Pressure 4-6 Bar
Protection Class IP-65

9 avcon


# Pneumatic Angle Seat Valve 6812M/6212M Series

## **6812M-Casted SPG, Pneumatically Operated Angle Seated Valve**




Valve Model No	Siz	е	Actuator (normally	Flow Kv	Pressure	Weight				Dimen	sions			
PTFE	INCH	ММ	close)	M₃Hr	Kg/cm2	Kgs	Α	С	D	Н	L1	L2	FTF	THK
6812M166/T/SPG	1/2"	15	K-51	4.3	12.0	0.97	182	157	44	12.7	72.5	32.5	89	1.65
6812M216/T/SPG	3/4"	20	K-51	5.4	12.0	1.18	187	162	44	19.05	90.5	29.5	120	1.65
6812M266/T/SPG	1"	25	K-51	13.6	10.0	1.83	209	180	44	25.4	96	29	150	1.65
6812M416/T/SPG	1 1/2"	40	K-81	31.8	10.0	2.35	265	226	60.5	38.1	133	27	160	1.65
6812M516/T/SPG	2"	50	K-81	54.8	10.0	3.27	270	226	60.5	50.8	148	32	180	1.65
6812M666/T/SPG	2 1/2"	65	K-121	59	6	5.8	347	295	87	63.5	240	50	290	1.65
6812M816/T/SPG	3"	80	K-121	65	5	6.0	375	320	87	76.2	255	55	310	1.65

## **6812M-Casted TC, Pneumatically Operated Angle Seated Valve**




	Valve Model No	Siz	:e	Actuator (normally	Flow Kv	Pressure	Weight			Dir	nensio	ns		
_	PTFE	INCH	MM	close)	M₃Hr	Kg/cm2	Kgs	Α	С	D	L1	L2	FTF	TD
	6812M167/T/TC	1/2"	15	K-51	4.3	12.0	0.97	182	157	44	71.6	30	101.6	25
	6812M217/T/TC	3/4"	20	K-51	5.4	12.0	1.18	187	162	44	71.6	30	101.6	50.4
	6812M267/T/TC	1"	25	K-51	13.6	10.0	1.83	209	180	44	81.3	33	114.3	50.4
	6812M417/T/TC	1 1/2"	40	K-81	31.8	10.0	2.35	270	226	60.5	112.7	27	139.7	50.4
	6812M517/T/TC	2"	50	K-81	54.8	10.0	3.27	277	230	60.5	126.7	32	158.7	63.9
	6812M667/T/TC	2 1/2"	65	K-121	59	6	5.8	360	290	87	143.8	50	193.8	77.4
	6812M817/T/TC	3"	80	K-121	65	5	6.0	375	320	87	167.3	55	222.3	90.9

## **6212M - Casted SPG, Hand Operated Angle Seated Valve**



ı	Valve Model No	Siz	e	Flow Kv	Pressure	Weight			Dim	ensio	ns
	PTFE	INCH	MM	M₃Hr	Kg/cm2	Kgs	Α	ØB	С	D	TD
	6212M166/T/SPG	1/2"	15	4.3	12.0	0.97	89	12.7	9.4	105	25
	6212M216/T/SPG	3/4"	20	5.4	12.0	1.18	120	19.05	15.75	122	50.4
	6212M266/T/SPG	1"	25	13.6	10.0	1.83	150	25.4	22.1	140	50.4
	6212M416/T/SPG	1 1/2"	40	31.8	10.0	2.35	160	35.1	34.8	140	50.4
	6212M516/T/SPG	2"	50	54.8	10.0	3.27	180	50.8	47.5	164	63.9

## **6212M - Casted TC, Hand Operated Angle Seated Valve**



Valve Model No	Size		Flow Kv	Pressure	Dimensions					
PTFE	INCH	MM	M₃Hr	Kg/cm2	Kgs	Α	ØB	ØC	Н	
6212M167/T/TC	1/2"	15	4.3	12.0	0.97	101.6	25	9.4	188	
6212M217/T/TC	3/4"	20	5.4	12.0	1.18	101.6	50.4	15.75	224	
6212M267/T/TC	1"	25	13.6	10.0	1.83	114.3	50.4	22.1	242	
6212M417/T/TC	1 1/2"	40	31.8	10.0	2.35	139.7	50.4	34.8	252	
6212M517/T/TC	2"	50	54.8	10.0	3.27	158.7	63.9	47.5	327	

10 avcou

## **Zero Dead Leg Block Valve**

#### **Design & Manufacturing Standards**

Design Standard AS
Testing Standard EN
Certification Standard EN

Surface Finish
Diaphragm Compliance

ASME BPE EN12266-1

EN-10204: 2.1,2.2,3.1 ASME BPE (SF4) USP CLASS VI Biological Reactivity Test In Vitro and In Vivo US FDA CFR177.1550(PTFE) US FDA CFR177.2600(EPDM)

#### **Benefits**

As per Customer's specific P&ID design

Compact design

Combination of many different Pipe diameters

Optimized drain ability Minimized dead leg

Reduced surface contact, hold up volume & cross

Contamination of the product

Reduction of multiple fittings, tubing & Piping in the

system

Reduced Qualification & Validation Documentation

#### **Applications**

Isolating Pipelines

Distribution

Point of Use

Sampling

Diverting, Mixing

Bypass Draining

Process Sterilization (SIP/CIP).

Dosing Feeding

#### **Operating Conditions**

Service Purified Water/ WIFI/ AIR/ GAS/ CIP/

SIP Intermediate and final products in the pharmaceutical & biotechnology

industries.

Operating Temp Ambient/85°C(EPDM) 155°C(PTFE)

Design Temp 150 ° C Operating Pr. 10 Kg/cm²(g) Design Pressure 12 Kg/cm²(g)/ FV

Vapour Yes

Viscosity cP As of service

Specific Gravity ~1

pH As of service

#### **Material of Construction**

Body ASTM A 182 SS 316L

Diaphragm Material EPDM or PTFE with EPDM backup

Studs/Nuts SS 304

#### **Pneumatic Actuator**

Type Linear

Construction Piston actuated (Single acting or double acting)

Material of construction Plastic-PA 6GF30 & SS 304 Hand Wheel Plastic-PA 6GF30 & SS 304

#### **Input for Pneumatic Actuator**

Air Signal Supply Bar (g) 4 to 6 BAR
Pneumatic Connections 1/4" BSP (F)
Fail Safe Position Normally Closed

Local Position Indicator Impact Resistant Thermoplastic Dome

#### **End Connections**

Inlet Spigot end as per BS 4825/ASME BPE Outlet Spigot end as per BS 4825/ASME BPE



#### **Finish**

Internal Electro polished to < 0.38 µm

Ra as per ASME BPE (SF4)

External Mirror Polished to 180 Grit

#### **Performance Characteristics**

Hydro Test Pressure Air Test Pressure Body 15.0 Kg/cm<sup>2</sup>(g) Seat 6Kg/cm<sup>2</sup>(g)

### Accessories

#### Solenoid Valves (optional)

Type 3/2 way, Normally close energies to open Coil voltage 230AC/24 DC/suitable for safe area

Air Supply Bar (g) 4 to 8

Material PA6/Brass/SS-304
Air Connection 1/4" BSPT (F)
Cable Entry Plug in connector PG-9

Area Classification Safe area, Weatherproof suitable for min. IP 65

#### Proximity Sensor (optional)

Proximity Switch

Make P&For equivalent

Type Inductive Proximity switches, 2 wire

Area classification Safe area, Weatherproof suitable for min. IP 65 No. of Switches 2 Nos. (1 No. Each for open & close positions)

Cable Entry 1/2" NPT (F) - 2 Nos. (1 No. Plugged)

#### Stroke Limiter (optional)

Function Actuators of diaphragm Valves

#### **Smart Positioner** (optional)

MOC PA6-GF30
Signal Input 4-20mA
Power supply 24V DC±10%
Ambient temperature 0-70°C
Supply Air pressure 4-6 Bar
Protection class IP-66
Explosion-proof class Ex nA II Ct4

11 avcon<sup>6</sup>

## **Multi Port Zero Dead Leg Block Valve**

## 4/2 Way, "T" Port Integrated Block Valve



## 4/2 Way, "L" Port Integrated Block Valve



## "L" Port ZDL Manual Block Valve



## "T" Port ZDL Manual Block Valve





#### MULTI PORT BLOCK VALVES

The technological evaluation in the Pharmaceuticals & Biotech industries, demands for stringent quality, Human safety & High efficiency from process system components.

To meet these growing demands, Avcon controls are offering engineering solutions to process piping by developing Customised multiport block valves. Our block valves provides minimum deadlegs, maximum drainability & simplify system Validation.

Our highly experience Design & development team is using highend design Software to assist in 3D modelling, fluid dynamics & finite element analysis **Advantages & Key Features.** 

- Block bodies are machined from forged blocks upto 2". Above 2" they are machined from the plates as per ASTM A 262/262M & ASTM A 240/240M.
- Material SS-316L as per ASME BPE.
- Customised multiport configurations for piping systems
- Compact design provides space optimisation
- Reduce hold up volume
- Simplified Validation by single block body
- Integrated design provides easy installation shorter CIP, enhanced cleanability, minimise system welds.
- Available with Avcon's standard diaphragms actuators

#### **Multi-port valve Application**

- Mixing
- Diverting
- Sampling
- Dumping
- Auto Switching
- CIP/SIP/Product Transfer



.2 **avcon**°

#### **Design & Manufacturing Standards**

Design Standard Testing Standard Certification Standard Surface Finish Diaphragm Compliance ASME BPE EN12266-1

EN-10204:2.1,2.2,3.1 ASME BPE (SF4) USP CLASS VI Biological Reactivity Test In Vitro and In

Vivo US FDA CFR177.1550(PTFE) US FDA CFR177.2600(EPDM)

#### **Benefits**

Very compact design, ideal for applications

where space is at a premium

Minimal deadlegs and optimized draining capabilities The valve body is machined from a single piece of block

Material (mono block-no welded components)

CIP/SIP capable and sterilisable

The internal surface finish of the valve body is available

Mechanically and/or electro polished is Ra 0.38  $\mu m$  Optimised flow geometry

Welding into the tank bottom is simplified by a welding neck (standard 6 mm)

Available with manual, pneumatic actuator

Optical and/or electrical position indicators are available as actuator instrumentation

#### **Applications**

Manual or air actuated outlet valve for sterile

vessel

. Bioreactor/Fermenter:vessels

Perfused bioreactors: retentate or permeate

vessels

Downstream: Outlet valve for process vessels in

the following production

Areas: separation, filtration, chromatography,

formulation

### **Operating Conditions**

Service Manual or air actuated outlet

valve for sterile vesel.

Operating Temp Ambient/85°C(EPDM) 155°C(PTFE)

Design Temp 150 ° C Operating Pr. 10 Kg/cm²(g) Design Pressure 12 Kg/cm²(g)/ FV

Vapour Yes

Viscosity cP As of service

Specific Gravity ~1

pH As of service

#### **Material of Construction**

Body ASTM A 182 SS 316L

Diaphragm Material EPDM or PTFE with EPDM backup

Studs/Nuts SS 304

#### **Pneumatic Actuator**

Type Linear

Construction Piston actuated (Single

acting or double acting) Manual Plastic-PA 6GF30 & SS 304

Material of construction Plastic-PA 6GF30 & SS 304 Hand Wheel Plastic-PA 6GF30 & SS 304

#### **Input for Pneumatic Actuator**

Air Signal Supply Bar (g) 4 to 6 BAR
Pneumatic Connections 1/4" BSP (F)
Fail Safe Position Normally Closed

Local Position Indicator Impact Resistant Thermoplastic Dome

#### **End Connections**

Inlet Spigot end as per BS 4825/ASME BPE Outlet Spigot end as per BS 4825/ASME BPE

## Flush Bottom Valve 6829 FBV (Radial) Tank Sampling Valve 6226 TSV (Radial)



#### **Finish**

Internal Electro polished to < 0.38 µm

Ra as per ASME BPE (SF4)

External Mirror Polished to 180 Grit

#### **Performance Characteristics**

Hydro Test Pressure Body 15.0 Kg/cm²(g)
Air Test Pressure 6Kg/cm²(g)

Accessories

#### Solenoid Valves (optional)

Type 3/2 way, Normally close energies to open Coil voltage 230AC/24 DC/suitable for safe area

Air Supply Bar (g) 4 to 8

Material PA6/Brass/SS-304 Air Connection 1/4" BSPT (F) Cable Entry Plug in connector PG-9

Area Classification Safe area, Weatherproof suitable for min. IP 65

#### **Proximity Sensor** (optional)

Proximity Switch

Make P&For equivalent

Type Inductive Proximity switches, 2 wire

Area classification Safe area, Weatherproof suitable for min. IP 65 No. of Switches 2 Nos. (1 No. Each for open & close positions)

Cable Entry 1/2" NPT (F) - 2 Nos. (1 No. Plugged)


#### Stroke Limiter (optional)

Function Actuators of diaphragm Valves

#### **Smart Positioner** (optional)

MOC PA6-GF30
Signal Input 4-20mA
Power supply 24V DC±10%
Ambient temperature 0-70°C
Supply Air pressure 4-6 Bar
Protection class IP-66
Explosion-proof class Ex nA II Ct4

13 avcor



#### TANK BOTTOM VALVES

The AVCON Tank Bottom Valves are designed for applications in the aseptic process industry offering a pocket-free interior surface, minimized sump, eliminating entrapment areas and minimized flow resistance, hence eliminating potential process contamination.

The tank valve body is machined from solid bar stock material 316L as per ASME BPE. The standard design offers one valve port outlet. Different options are available for sampling, sterilization and multi-outlet configurations for the customized solutions.

To minimizes the hold up volume we recommend to weld TBV directly to vessel. Tank Bottom valves are typically used for tank discharge, draining, sampling, cleaning and / or sterilizing, rinsing and isolation of down stream processing.

#### **Features:**

- Valve body machined from a solid bar stock material.
- Material 316L as per ASME BPE.
- Minimized dead leg and internal
- Suitable for mounting with AVCON K-Series Pneumatic Actuator.

## **Accessories**

- Stainless Steal Actuator
- Feedback Position
- Positioner
- Spare Diaphragm

3) Manual TBV OR 1x Main Valve 1x Sample Valve

Outlet left

side right

Manual TBV 1x Main Valve 1x Sample Valve **Outlet Right** 

side left

- 1x CIP/SIP cleaning outlet valve left
- 5)Manual TBV
- 1x Main Valve
- 2x Outlet port for loop installation or as two access ports.

#### Sanitary Pneumatic SS Actuated Ball Valve - 2390(V) / 2390(R) Ball Valve - 2370(V) / 2370(R) **Butterfly Valve - 4322M(V) / 4322M(R)**

#### **Design & Manufacturing Standards**

Design Standard Testing Standard EN12266-1 Surface Finish ASME BPE (SF:2.4-1)

Conformity Standard EN:10204

#### Pneumatic Actuator

Linear Vertical Actuator / Rotary Actuator Type Actuator Body SS-304 SS-304 Bracket Signal Pressure 3.5 to 6 bar Pneumatic Connection 1/4" BSP

> **Ball Valve** 2370(V)

**Ball Valve** 2370M(R)



#### **Process Parameters**

Service Air/Water/Steam/Slightly Pneumatic Connection 1/4" BSP Working Temperature PTFE upto 180°C

Working Pressure 10 Bar

#### **Benefits**

Maintenance-free pneumatic Rotary & Linear actuator With "Normally Closed",

"Normally Open" or "Double Acting" Orientations. Available with the control unit & interface unit

Standard ISO pattern mounting pad Stainless steel Construction: Forged body & ball.

Two & Three Piece design with floating seals Optionally available with cavity-filled seat

FDA approved PTFE seat

High temperature & pressure rating

Low maintenance

Low torque

High-grade surface finish

Suitable for vacuum applications

CIP/SIP capable

#### **Applications**

Food & Drinks:

Soups & Sauces, Flavorings & Ingredients, Brewery, Wine,

Fats & Oils, Animal Oils.

Personal Care: Extracts, Face Creams & Lotions Perfumes Soaps. Hair Styling Gels & Liquids, Dyes

Chemical: Paints, Adhesives, Coatings, Oils & Lubricants, Detergents Emulsions, Fuels.

Utilities: Steam, CIP Flush, High Purity Water, Water

for Injection (WFI).

#### Valve Conditions

Valve Body SS.304/316/316L as per ASTM

A 276A/182A

**Body Pressure Rating** 16 Bars Surface Finish Internal: 0.4 Ra Max

**End Connections** Spigot & TC as per ASME BPE/

BS 4825

**Butterfly Valve** 4322(V)

**Butterfly Valve** 4322M(R)





#### **Process Parameters**

Air/Gas/Liquid Service

Working Temperature For EPDM: (-) 10°C to 90°C and

Silicon upto 120°C

Working Pressure 10 Bar

**Elastomer** 

EPDM/Silicon Seat Compliance

#### **Benefits**

Maintenance-free pneumatic Rotary & Linear actuator With "Normally Closed",

"Normally Open" or "Double Acting" orientations.

Available with the identical control unit, controls interface Standard ISO pattern mounting pad.

Construction: Forged SS-304, SS-316 or 316L body & disc.

FDA-approved seat materials: EPDM, Silicone. Sanitary surface finish, with Ra < 0.4 µm or EP (SF4)

Suitable for low and medium viscosity fluids.

100% hydraulic tested / 100% inspected.

#### **Applications**

CIP Lines

Food & Beverage

Brewing, Distilling, and Wine Making

Dairy


Pharmaceuticals

Chemical Processing


Flavor and Fragrance

15 avcon

## **Sanitary Type Pneumatic Vertical Actuated Ball Valve - 2370M(V)**





CAT. NO.	SIZE	ID	PIPE OD	THK	Q	Н	Р	SQ.	L
2370M166	1/2"	9.4	12.7	1.65	55	53.25	275.25	11	95
2370M216	3/4"	15.7	19.05	1.65	60	54.5	276.5	11	106
2370M266	1"	22.1	25.4	1.65	75	65	287	11	118
2370M416	1-1/2"	34.8	38.1	1.65	100	72	294	11	145
2370M516	2"	47.5	50.8	1.65	120	82.5	304.5	14	165
2370M666	2-1/2"	63	66.3	1.65	140	99	321	17	191




CAT. NO.	SIZE	PIPE OD	TC OD	Н	Р	SQ.	L
2370M167	1/2"	9.3	25	53.25	275.25	11	69.2
2370M217	3/4"	15.6	35	54.5	276.5	11	80.6
2370M267	1"	21.7	50.4	65	287	11	93
2370M417	1-1/2"	34	50.4	72	294	11	105
2370M517	2"	47.5	63.9	82.5	304.5	14	114
2370M667	2-1/2"	63	77.4	99	321	17	144

## Sanitary Type Pneumatic Vertical Actuated Butterfly Valve - 4322(V)









CAT. NO.	SIZE	PIPE ID	TC OD	В	С	D	Е	F	G	SQ
4322M267	1"	22.1	50.4	75	22	60	86	123.5	343.5	11
4322M417	1 1/2"	34.8	50.4	77	22	60	87.5	126	348	11
4322M517	2"	47.5	63.9	85	22	60	91.5	134	356	11
4322M667	2 1/2"	60.2	77.4	108	22	60	104	158	380	11
4322M817	3"	72.9	90.9	122	22	72	114	178	400	14
4322M1017	4"	97.3	118.9	153	22	80	127.5	205	427	14
4322M1517	6"	146.8	166.9	200	25	85	152	252	474	14
4322M2017	8"	197.6	217.5	250	28	88	177	302	524	17


16 avcou


## Sanitary Type Pneumatic Rotary Actuated Ball Valve - 2370M (R)






CAT. NO.	ACTUATOR MODEL	SIZE	PIPE ID	PIPE OD	THK	Q	Н	Р	SQ.	L	K	J
2370M166	ACT-063SR06	1/2"	9.4	12.7	1.65	55	53.25	175	11	69.2	173	82
2370M216	ACT-075SR06	3/4"	15.7	19.05	1.65	60	54.5	185	11	80.6	221	97
2370M266	ACT-090SR06	1"	22.1	25.4	1.65	75	65	210.5	11	93	260	104
2370M416	ACT-090SR07	1-1/2"	34.8	38.1	1.65	100	72	222.5	11	105	260	104
2370M516	ACT-115SR06	2"	47.5	50.8	1.65	120	82.5	241	14	114	312	132
2370M666	ACT-125SR06	2-1/2"	63	66.3	1.65	140	99	270	17	144	368	149






CAT. NO.	ACTUATOR MODEL	SIZE	PIPE ID	TC OD	L	K	J	Н	Р	SQ.
2370M167	ACT-063SR06	1/2"	9.3	25	69.2	173	82	53.25	175	11
2370M217	ACT-075SR06	3/4"	15.6	35	80.6	221	97	54.5	185	11
2370M267	ACT-090SR06	1"	21.7	50.4	93	260	104	65	210.5	11
2370M417	ACT-090SR07	1-1/2"	34	50.4	105	260	104	72	222.5	11
2370M517	ACT-115SR06	2"	47.5	63.9	114	312	132	82.5	241	14
2370M667	ACT-125SR06	2-1/2"	63	77.4	144	368	149	99	270	17

## Sanitary Type Pneumatic Rotary Actuated Butterfly Valve - 4322M (R)



CAT. NO.	ACTUATOR MODEL	SIZE	PIPE ID	PIPE OD	А	В	С	D	Н	J	G	SQ
4322M416	ACT-050D	1-1/2"	34.8	38.1	1.65	77	22	60	157	70	212.5	11
4322M516	ACT-063D	2"	47.5	50.8	1.65	85	22	60	173	82	239	11
4322M666	ACT-063D	2-1/2"	60.2	63.5	1.65	108	22	60	173	82	263	11
4322M816	ACT-063D	3"	72.9	76.2	1.65	122	22	72	173	82	283	14
4322M1016	ACT-090D	4"	97.3	101.6	2.11	153	22	80	260	104	340	14
4322M1516	ACT-0100D	6"	146.8	152.4	2.77	200	25	85	290	122	399	14
4322M2016	ACT - 0115D	8"	197.6	203.2	2.77	250	28	88	312	132	472	17



CAT. NO.	ACTUATOR MODEL	SIZE	PIPE ID	TC OD	В	С	D	Н	J	G	SQ
4322M417	ACT-050D	1-1/2"	34.8	50.4	77	22	60	157	70	212.5	11
4322M517	ACT-063D	2"	47.5	63.9	85	22	60	173	82	239	11
4322M667	ACT-063D	2-1/2"	60.2	77.4	108	22	60	173	82	263	11
4322M817	ACT-063D	3"	72.9	90.9	122	22	72	173	82	283	14
4322M1017	ACT-090D	4"	97.3	118.9	153	22	80	260	104	340	14
4322M1517	ACT-0100D	6"	146.8	166.9	200	25	85	290	122	399	14
4322M2017	ACT - 0115D	8"	197.6	217.5	250	28	88	312	132	472	17

17 avcon

is highly motivated team of Technical, Managerial Professionals with strong commitment to excellence in Manufacturing / Supply of Automated Valves, which includes: SOLENOID MOTORISED PNEUMATICALLY OPERATED VALVES, used in diversed field where automation is needed.

## Motto is "Quality + Technology + Performance."

This ensures "The worth of money you spend in buying Automated Valves." Pioneer spirit at AVCON to design & develop Industrial duty Automated Valves, resulted in Manufacturing most comprehensive complete range of Automated Valves. Recognized everywhere and yet economically priced.







Office & Works Unit -II







Office & Works Unit -IV

## **Diaphragm Valves**

For Aseptic & Hygienic Application



## Our other most popular range of products



